

Certificate of Analysis

Product Name	Proteose peptone	
Cat Number	P8960	
Lot Number	2507050002	
Storage	RT	
Test	Specification	Result
Appearance	White to light yellow powder	Light yellow powder
Total nitrogen	≥ 12.0%	13.10%
Amino nitrogen	1.5-4.5%	3.60%
Loss on drying	≤ 6.0%	3.69%
Residue on ignition	≤ 9.0%	8.09%
Test Date	05 Jul 2025	
Recommended Retest Period	5 Years	
The results meet the quality requirements		
QC Checker	Niu Hao Ran	QC Manager



Quan Jing Yu

