

Certificate of Analysis

Product Name	Ficoll Plus 1.090		
Cat Number	P4370		
Lot Number	2500040005		
Storage	RT◆		
Test	Specification	Result	
Appearance	Colorless transparent liquid	Colorless transparent liquid	
pH	7.0-7.5	7.31	
Osmotic Pressure	290-350 mOsmol/kg	290 mOsmol/kg	
Density	1.090±0.001g/mL	Qualified	
Endotoxin	<5 EU/ml	<5 EU/ml	
Mycoplasma Detection	Negative	Negative	
Bacterial Detection	Sterile	Sterile	
Test Date	29 Apr 2025		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan Jing Yu

