

Certificate of Analysis

| | | |
|---|--|------------|
| Product Name | Glycolic acid | |
| CAS. | 79-14-1 | |
| Cat Number | G8841 | |
| Lot Number | 1550522001 | |
| Molecular Formula | C ₂ H ₄ O ₃ | |
| Molecular Weight | 76.05 | |
| Storage | RT | |
| Test | Specification | Result |
| Appearance | White to off-white solid | Qualified |
| Purity | ≥99% | 99% |
| Solubility (100mg/ml, H ₂ O) | Colorless transparent solution | Qualified |
| Test Date | 05 Dec 2025 | |
| Recommended Retest Period | 3 Years | |
| The results meet the quality requirements | | |
| QC Checker | Liu Juan | QC Manager |

