

Certificate of Analysis

Product Name	Trifluridine		
CAS.	70-00-8		
Cat Number	T8271		
Lot Number	504C032		
Molecular Formula	C ₁₀ H ₁₁ F ₃ N ₂ O ₅		
Molecular Wt	296.20		
Storage	2-8℃		
Test	Specification	Result	
Appearance	White to off-white crystalline powder	Off-white crystalline powder	
Purity	>98%	Qualified	
Solubility (5mg/100ul, DMSO)	Colorless transparent solution	Colorless transparent solution	
Loss on drying	≤1.0%	0.51%	
Test Date	19 Oct 2024		
Recommended Retest Period	3 Years		
The results meet the quality requirements			
QC Checker	Niu Hao Ran	QC Manager	Quan Jing Yu

