

Certificate of Analysis

Product Name	Dialysis Membranes(8000-14000D)		
Cat Number	YA1074		
Lot Number	24240702003		
Storage	RT		
Test	Specification	Result	
Tightness	No leakage and good seal	Qualified	
Temperature resistance	Can be boiled and high pressure	Qualified	
Dialysis experiment	Qualified	Qualified	
Materials	Regenerated cellulose	Qualified	
Test Date	27 Nov 2024		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan ling Yu

