

Certificate of Analysis

Product Name	Succinic anhydride	
CAS.	108-30-5	
Cat Number	S5620	
Lot Number	2512080003	
Molecular Formula	C ₄ H ₄ O ₃	
Molecular Weight	100.07	
Storage	2-8°C	
Test	Specification	Result
Appearance	White solid	Qualified
Purity	≥98%	98%
Solubility(100mg/ml ,DMSO)	Colorless clear solution	Qualified
Test Date	08 Dec 2025	
Recommended Retest Period	2 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager Wei Yan Hui

