



TIANJIN FCOMBIO SCIENTIFIC CO., LTD

Address: No.21 Zaolin Avenue, Xiqing District, Tianjin, China
Tel: 0086 13820252415 Web: www.fcombio.com Email: penny@fcombio.com

Product Quality Inspection Report

Product Name	Centrifuge tube	Part NO.	CT005-01S	Lot No.	2025032002
MF Date	2025.03.20	EXP Date	2030.03.19		
No.	Test item	Test content/inspection basis			Test result
1	Appearance	Black spots, bubbles, pull marks, etc.			meet the requirements
2	Dimensions	Meet the specifications of the drawings			meet the requirements
3	Printing	The printed content is correct and clear, without ink shortage or ghosting Printing is firm, 3M tape is pasted 3 times without deinking			NA
4	Hydrophilic	After hydrophilic TC treatment, measure the contact angle			NA
5	Centrifugal force	Comply with SIP regulation centrifugal force value			meet the requirements
6	Tightness test	Vacuum -30kpa~-29kpa, no leakage in 3 minutes			meet the requirements
7	Tightness test	Vacuum -30kpa~-29kpa, no leakage in 3 minutes			meet the requirements
8	Low retention test	10% Tween solution Low retention residue test			meet the requirements
9	Sterilization	The irradiation label is red, with "Certificate of Irradiation"			meet the requirements
10	Endotoxin testing	"Bacterial Endotoxin Inspection Specification"			meet the requirements
11	Cell culture testing	"Specifications for Performance Inspection of Cell Culture Products"			NA
12	Rnase detection	"Nuclease Test Specification"			meet the requirements
13	Dnase detection	"Nuclease Test Specification"			meet the requirements
Comprehensive judgment		OK			

Note: "inspection result" can be filled: meets the requirements, does not meet the requirements or NA

APPROVED BY :

徐骏

TITLE : QC MANAGER

DATE : 2025.03.25

