

Certificate of Analysis

Product Name	Cysteamine		
CAS.	60-23-1		
Cat Number	C7820		
Lot Number	2509170001		
Molecular Formula	C ₂ H ₇ NS		
Molecular Wt	77.15		
Storage	RT		
Test	Specification	Result	
Appearance	White crystalline powder	White crystalline powder	
Solubility (0.1g/ml, H ₂ O)	Colorless and transparent solution	Colorless and transparent solution	
Purity	Assay ≥95%	97%	
Melting Range	95-100℃	Qualified	
Test Date	17 Sep 2025		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan Jing Yu

