

Certificate of Analysis

Product Name	Diphenyliodonium Chloride	
CAS.	1483-72-3	
Cat Number	D9550	
Lot Number	2512060002	
Molecular Formula	C ₁₂ H ₁₀ ClI	
Molecular Weight	316.57	
Storage	RT	
Test	Specification	Result
Appearance	White to light yellow powder	White to light yellow powder
Assay	98.00%	Qualified
Solubility	Dissolved in methanol	Qualified
Test Date	06 Dec 2025	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager

