

Certificate of Analysis

Product Name	L-Glutamicacid sodiumsalt		
CAS.	6106-04-3		
Cat Number	G8010		
Lot Number	34230913007		
Molecular Formula	$C_5H_8NNaO_4 \cdot H_2O$		
Molecular Weight	187.13		
Storage	RT◆		
Test	Specification	Result	
Appearance	White crystalline or crystalline powder	Qualified	
Assay	≥98%	98.5%	
Solubility (40mg/ml, H ₂ O)	Colorless transparent solution	Colorless transparent solution	
Test Date	05 Nov 2024		
Recommended Retest Period	3 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan Jing Yu