

## Certificate of Analysis

Product Name	PMSF (100mM)		
CAS.	329-98-6		
Cat Number	P0100		
Lot Number	2500040004		
Molecular Formula	C <sub>7</sub> H <sub>7</sub> SO <sub>2</sub> F		
Storage	-20℃ ◆		
Test	Specification	Result	
Appearance	Colorless transparent liquid	Colorless transparent liquid	
Concentration	100mM	100mM	
Function Detection	Protease Inhibitor	Qualified	
Bacterial Detection	sterile	Qualified	
Protease	None	Qualified	
RNase	None	Qualified	
Test Date	25 Apr 2025		
Recommended Retest Period	1 Year		
The results meet the quality requirements			
QC Checker	Niu Hao Ran	QC Manager	Quan Jing Yu

