

Certificate of Analysis

Product Name	3, 5-Diaminobenzoic Acid	
CAS.	535-87-5	
Cat Number	D1050	
Lot Number	2509180002	
Molecular Formula	C ₇ H ₈ N ₂ O ₂	
Molecular Wt	152.15	
Storage	RT	
Test	Specification	Result
Appearance	Crystalline powder	Crystalline powder
Melting Range	235-238℃	235℃
Test Date	18 Sep 2025	
Recommended Retest Period	2 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager 质检章 Quan Jing Yu

