

Certificate of Analysis

Product Name	PMA		
CAS.	16561-29-8		
Cat Number	P6741		
Lot Number	2550308004		
Molecular Formula	C ₃₆ H ₅₆ O ₈		
Molecular Weight	616.8		
Storage	-20℃		
Test	Specification	Result	
Appearance	White powder	White powder	
Assay	≥98.0%	99.6%	
Solubility (10mg/ml, DMSO)	Colorless and transparent solution	Colorless and transparent solution	
Loss on Drying	≤2.0%	0.5%	
Test Date	05 Jul 2025		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan Jing Yu

