


Certificate of Analysis

Product Name	Dansyl Chloride	
CAS.	605-65-2	
Cat Number	D9580	
Lot Number	3540702001	
Molecular Formula	C ₁₂ H ₁₂ ClNO ₂ S	
Molecular Weight	269.75	
Storage	2-8℃	
Test	Specification	Result
Appearance	Yellow powder	Yellow powder
Assay	≥99.0%	Qualified
Solubility (20mg/ml, CHCl ₃)	Orange clear solution	Orange clear solution
Test Date	29 Apr 2025	
Term of validity	3 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager



Wei Yan Hui

