

Certificate of Analysis

Product Name	Lysostaphin	
CAS.	9011-93-2	
Cat Number	L9070	
Lot Number	24230816004	
Storage	-20℃	
Test	Specification	Result
Appearance	White or light yellow powder	Qualified
Purity (SDS-PAGE)	>80%	Qualified
Solubility (50mmol/L pH 7.5 Tris-HCl)	Colorless or light yellow transparent solution	Qualified
Enzyme Activity	>1000U/mg	2000U/mg
Test Date	19 Dec 2024	
Recommended Retest Period	2 Years	
The results meet the quality requirements		
QC Checker	Gao Xue	QC Manager Wei Yan Hui

