

Certificate of Analysis

Product Name	Dextran 5000		
CAS.	9004-54-0		
Cat Number	D7270		
Lot Number	2541012001		
Storage	RT		
Test	Specification	Result	
Appearance	White to off-white powder	Qualified	
Solubility (50mg/ml, H ₂ O)	Colorless clear solution	Colorless clear solution	
pH (1%, H ₂ O)	4.5-7.0	5.28	
Loss on drying	≤7.0%	4.96%	
Test Date	15 Sept 2025		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Niu Hao Ran	QC Manager	Quan Jing Yu

