

Certificate of Analysis

Product Name	HP-20 Macroporous Adsorption Resin	
Cat Number	M0043	
Lot Number	1551017003	
Storage	RT	
Test	Specification	Result
Appearance	White opaque particles	Qualified
Water Content	65%—75%	68.0%
Specific Surface Area	550-600	Qualified
Particle size	0.3-1.25	Qualified
Wet True Density(g/ml)	1.03-1.07	1.05
Wet Apparent Density(g/ml)	0.68-0.75	0.70
Test Date	20 Oct 2025	
Recommended Retest Period	1 Year	
The results meet the quality requirements		
QC Checker	Niu Hao Ran	QC Manager Quan Jing Yu

