

## Certificate of Analysis

Factory Name: Tianjin Fcombio Scientific Co.,Ltd  
Product Name: Syringe Filter (Non sterile) SF25TH-045  
Size: 25mm 0.45µm Hydrophilic PTFE  
Housing Material: PP (Polypropylene)  
Lot No.: 210805

Each lot has been sampled, tested and released by the Quality Assurance Department according to ISO9001:2016 for the following characteristics:

Specification	Acceptable limits	Test Data
Filtration Area	3.8 cm2	√
Holdup Volume (ul)	125 ul	√
Wettability	Hydrophilic	√
Volume Throughput (ml)	100	√
Inlet/Outlet	Female luer Lock/ Male Luer Slip	√
Maximum Operation Temperature(°C)	100	100
99.5% EtOH Bubble Point at 23 °C	0.85bar	√
Bursting Pressure	>10 bar	√
Cleanliness level	No spot, no Impurity	√
Housing Integrity Test	No leaking under pressure	√
Housing appearance	No Flashes, no deform	√

The product corresponded to the requirements specified by the technical condition upon quality test and permitted to leave the factory

Certificate Approved By:Wandan

Issue Date: 10th Mar 2025

