

Certificate of Analysis

Product Name	N-Hydroxysuccinimide NHS	
CAS.	6066-82-6	
Cat Number	H8220	
Lot Number	1550704004	
Molecular Formula	C ₄ H ₅ NO ₃	
Molecular Wt	115.09	
Storage	2-8℃	
Test	Specification	Result
Appearance	White crystalline powder	White crystalline powder
Assay	≥99.0%	99.4%
Solubility (100mg/ml, H ₂ O)	Colorless clear transparent solution	Colorless clear transparent solution
Water content	≤2.0%	0.33%
Test Date	10 Nov 2025	
Recommended Retest Period	1 Year	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager

索宝科技

有限公司

检测

合格

2025.11.10

Quan Jing Yu

