

## Certificate of Analysis

|                                           |                                 |                     |
|-------------------------------------------|---------------------------------|---------------------|
| Product Name                              | Pectin                          |                     |
| CAS.                                      | 9000-69-5                       |                     |
| Cat Number                                | P8030                           |                     |
| Lot Number                                | 2651124001                      |                     |
| Storage                                   | RT                              |                     |
|                                           |                                 |                     |
| Test                                      | Specification                   | Result              |
| Appearance                                | White to brownish yellow powder | Light yellow powder |
| Moisture Content                          | ≤10%                            | 7.9%                |
| Galacturonic Acid                         | ≥58.0%                          | 83%                 |
| Methoxy Group                             | ≥6%                             | 8.0%                |
| Test Date                                 | 22 Jun 2026                     |                     |
| Recommended Retest Period                 | 2 Years                         |                     |
| The results meet the quality requirements |                                 |                     |
| QC Checker                                | Liu Juan                        | QC Manager          |



北京宝康生物科技有限公司  
Quen Jing Yu

