

Certificate of Analysis

Product Name	Sephadex LH-20	
Cat Number	S8111	
Lot Number	1540801003	
Storage	RT	
Test	Specification	Result
Appearance	White particles	Qualified
Average particle size	100±10um	Qualified
Swelling Volume	4±0.2mL/g	Qualified
Moisture Content	<3%	2.7%
Test Date	15 Sep 2025	
Recommended Retest Period	1 Year	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager 质检章 Guanying Yu

