

## Certificate of Analysis

Product Name	Goat Anti-Rat IgG/Solar Fluor 350		
Cat Number	K1032G-AF350		
Lot Number	1225110001		
Host	Goat		
Storage	-20℃		
Test	Specification	Result	
Appearance	Liquid	Qualified	
Recommended dilution	IF=1:50-500	Qualified	
Purification method	Affinity purified by Protein A	Qualified	
Test Date	11 Dec 2025		
Recommended Retest Period	1 Year		
The results meet the quality requirements			
QC Checker	Niu Hao Ran	QC Manager	Quan Jing Yu

