

Certificate of Analysis

Product Name		PPLO Broth Base	
Cat Number		LA7080	
Lot Number		20241106	
Storage		RT	
Test		Specification	Result
Appearance		Yellow uniform powder	Yellow uniform powder
Solubility (H₂O)		Bright yellow clear solution	Bright yellow clear solution
pH		7.6-8.0	7.78
Moisture content		≤6.0%	3.9%
Microbiological Performance (35-37 °C ,7-14 d)	Test bacteria	Inoculum (cfu)	Test medium
	<i>Mycoplasma pneumonia</i> ATCC 15531	10 ⁻⁸	Growing well
	<i>Mycoplasma orale</i> ATCC23714	10 ⁻⁴	Growing well
Recommended Retest Period		3 Years	
The results meet the quality requirements			
QC Checker	Niu Hao Ran	QC Manager	Quan Jing Yu

