

Certificate of Analysis

Product Name	50×TAE buffer	
Cat Number	T1060	
Lot Number	2600030001	
Storage	RT	
Test	Specification	Result
Appearance	Clear colorless solution	Clear colorless solution
pH（1×, 25℃）	8.2-8.4	8.25
Deoxyribonuclease	None detected	None detected
DNA Electrophoresis	Qualified	Qualified
Test Date	23 Mar 2026	
Recommended Retest Period	2 Years	
The results meet the quality requirements		
QC Checker	Niu Hao Ran	QC Manager



Quen Jing Yu

