

Certificate of Analysis

Product Name	Monensin Sodium	
CAS.	22373-78-0	
Cat Number	M8670	
Lot Number	3641010002	
Molecular Formula	C ₃₆ H ₆₁ NaO ₁₁	
Molecular Weight	692.85	
Storage	RT	
Test	Specification	Result
Appearance	Off white powder	Qualified
Solubility（50mg/ml EtOH）	Colorless and transparent solution	Qualified
Purity	> 90.0%	Qualified
Test Date	08 May 2026	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager Quan Jing Yu

