

Certificate of Analysis

Product Name	Phytigel plantcell	
CAS.	71010-52-1	
Cat Number	P8170	
Lot Number	3541118003	
Storage	RT	
Test	Specification	Result
Appearance	White or off-white powder	White powder
Assay	85.0-108.0%	Qualified
pH	5-8	Qualified
Gel strength	≥850g/cm ²	1080g/cm ²
Transmittance	≥78%	87.5%
Loss on drying	≤15.0%	10.0%
Granularity (passing 60 mesh sieve)	≥92%	Qualified
Test Date	01 Jul 2025	
Recommended Retest Period	3 Years	
The results meet the quality requirements.		
QC Checker	Liu Juan	QC Manager



质检章

