

Certificate of Analysis

Product Name	DTE	
CAS.	6892-68-8	
Cat Number	D8421	
Lot Number	34240702001	
Molecular Formula	C ₄ H ₁₀ O ₂ S ₂	
Molecular Weight	154.25	
Storage	2-8℃	
Test	Specification	Result
Appearance	White powder	Qualified
Purity	> 99%	Qualified
Solubility (50mg/ml, H ₂ O)	Colorless and transparent solution	Qualified
Test Date	03 Jul 2024	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Niu Hao Ran	QC Manager



