

Certificate of Analysis

Product Name	Carbamazepine		
CAS.	298-46-4		
Cat Number	C5920		
Lot Number	2507070501		
Molecular Formula	C ₁₅ H ₁₂ N ₂ O		
Molecular Weight	236.27		
Storage	2-8℃		
Test	Specification	Result	
Appearance	White or almost white crystalline powder	Qualified	
Assay	≥98%	98%	
Solubility (10mg/ml, Ethanol)	Colorless and transparent solution	Colorless and transparent solution	
Test Date	07 Jul 2025		
Recommended Retest Period	2 Years		
The results meet the quality requirements			
QC Checker	Liu Juan	QC Manager	Quan Jing Yu

