

Certificate of Analysis

Product Name	Anti-CSDE1 Polyclonal Antibody	
Cat Number	K109309P	
Lot Number	8123202005	
Host	Rabbit	
Subtype	IgG	
Storage	-20°C	
Test	Specification	Result
Appearance	Liquid	Liquid
Recommended dilution	WB 1:8000-12000 IHC 1:50-200 IF 1:50-200.	Qualified
Test Date	09 Aug 2023	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Niu Hao Ran	QC Manager Quan Jing Yu

