

Certificate of Analysis

Product Name	Fluconazole	
CAS.	86386-73-4	
Cat Number	F8380	
Lot Number	3530818002	
Molecular Formula	C ₁₃ H ₁₂ F ₂ N ₆ O	
Molecular Weight	306.27	
Storage	RT	
Test	Specification	Result
Appearance	White powder	White powder
Purity	≥98%	98%
Solubility (50mg/ml, MeOH)	Colorless transparent solution	Colorless transparent solution
Test Date	06 Aug 2025	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager Quan Jing Yu

