

## Certificate of Analysis

<b>Product Name</b>	Casein Enzymatic Hydrolysate	
<b>CAS.</b>	91079-40-2	
<b>Cat Number</b>	C8210	
<b>Lot Number</b>	4550905006	
<b>Storage</b>	RT	
<b>Test</b>	<b>Specification</b>	<b>Result</b>
<b>Appearance</b>	Yellow to light brown powder	Qualified
<b>Solubility (20mg/ml, H<sub>2</sub>O)</b>	Transparent or slightly cloudy yellow liquid	Transparent yellow liquid
<b>Loss on Drying</b>	≤5.0%	3.4%
<b>Nitrogen</b>	≥12.0%	13.0%
<b>Amino Content</b>	6.0-7.0%	6.5%
<b>Test Date</b>	12 Sep 2025	
<b>Recommended Retest Period</b>	4 Years	
The results meet the quality requirements		
<b>QC Checker</b>	Liu Juan	<b>QC Manager</b> Quan Jing Yu

