

Certificate of Analysis

Product Name	Sodium Lauroyl sarcosine	
CAS.	137-16-6	
Cat Number	S8100	
Lot Number	3550211002	
Molecular Formula	C ₁₅ H ₂₈ NNaO ₃	
Molecular Wt	293.39	
Storage	RT	
Test	Specification	Result
Appearance	White powder	White powder
Purity	≥95.0%	Qualified
Solubility (10%, H ₂ O)	Colorless clear solution	Colorless clear solution
pH (10%, H ₂ O)	7.5-8.5	Qualified
Sodium Laurate	≤4.0%	2.0%
Test Date	12 Feb 2025	
Recommended Retest Period	3 Years	
The results meet the quality requirements		
QC Checker	Liu Juan	QC Manager



Quan Jing Yu

